

ARCHIVE COPY

Naval Ocean Systems Center
San Diego, CA 92152-5000

Technical Document 1633
August 1989

Signal Processing Device to Control Microwave Output

San Diego State University Foundation

The views and conclusions contained in this report are those of the contractors and should not be interpreted as representing the official policies, either expressed or implied, of the Naval Ocean Systems Center or the U.S. Government.

Approved for public release; distribution is unlimited.

ARCHIVE COPY

NAVAL OCEAN SYSTEMS CENTER

San Diego, California 92152-5000

E. G. SCHWEIZER, CAPT, USN
Commander

R. M. HILLYER
Technical Director

ADMINISTRATIVE INFORMATION

This work was performed for the Naval Medical Research and Development Command, National Naval Medical Center, Bethesda, MD 20814-5044, under program element number 0603706N. Contract N66001-85-D-0203 was carried out by San Diego State University Foundation, 5178 College Avenue, San Diego, CA 92182-1900, under the technical coordination of E. Lindner, Chemistry/Biochemistry Branch, Code 521, Naval Ocean Systems Center.

Released by
E. Lindner, Head
Chemistry/Biochemistry Branch

Under authority of
S. Yamamoto, Head
Environmental Sciences
Division

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE

REPORT DOCUMENTATION PAGE

1a REPORT SECURITY CLASSIFICATION UNCLASSIFIED			1b RESTRICTIVE MARKINGS	
2a SECURITY CLASSIFICATION AUTHORITY			3 DISTRIBUTION AVAILABILITY OF REPORT	
2b DECLASSIFICATION/DOWNGRADING SCHEDULE			Approved for public release, distribution is unlimited	
4 PERFORMING ORGANIZATION REPORT NUMBER(S)			5 MONITORING ORGANIZATION REPORT NUMBER(S) NOSC TD 1633	
6a NAME OF PERFORMING ORGANIZATION San Diego State University Foundation		6b OFFICE SYMBOL (if applicable)	7a NAME OF MONITORING ORGANIZATION Naval Ocean Systems Center	
6c ADDRESS (City, State and ZIP Code) 5178 College Avenue San Diego, California 92182-1900		7b ADDRESS (City, State and ZIP Code) San Diego, California 92152-5000		
8a NAME OF FUNDING SPONSORING ORGANIZATION Naval Medical Research & Development Center		8b OFFICE SYMBOL (if applicable) NMRDC	9 PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER N66001-85-D-0203	
8c ADDRESS (City, State and ZIP Code) National Naval Medical Center Bethesda, MD 20814-5044		10 SOURCE OF FUNDING NUMBERS		
		PROGRAM ELEMENT NO 0603706N	PROJECT NO FA24	TASK NO 560-EE90
		AGENCY ACCESSION NO DN787 507		
11 TITLE (Include Security Classification) SIGNAL PROCESSING DEVICE TO CONTROL MICROWAVE OUTPUT				
12 PERSONAL AUTHOR(S) J. G. Pinto, Ph D				
13a TYPE OF REPORT Final	13b TIME COVERED FROM TO	14 DATE OF REPORT (Year, Month, Day) August 1989	15 PAGE COUNT 11	
16 SUPPLEMENTARY NOTATION				
17 COSATI CODES			18 SUBJECT TERMS (Continue on reverse if necessary and identify by block number)	
FIELD	GROUP	SUB-GROUP	microwave oven signal processing device	
19 ABSTRACT (Continue on reverse if necessary and identify by block number) The development of an electronic device to control the operation of a commercial microwave oven is discussed. This device, when installed in conjunction with the existing circuitry of SHARP MICROWAVE OVEN (model R-9524), is capable of automatically advancing through a sequence of thawing recipes programmed and stored in the memory bank of the oven. The device therefore eliminates or minimizes human operator action needed in previous prototype versions of a blood-thawing device.				
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT <input type="checkbox"/> UNCLASSIFIED/UNLIMITED <input checked="" type="checkbox"/> SAME AS RPT <input type="checkbox"/> DTIC USERS			21 ABSTRACT SECURITY CLASSIFICATION UNCLASSIFIED	
22a. NAME OF RESPONSIBLE PERSON E. Lindner		22b. TELEPHONE (Include Area Code) (619) 553-2795	22c OFFICE SYMBOL Code 521	

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

TITLE A Signal Processing Device to Control Microwave Output.

ABSTRACT

The development of an electronic device to control the operation of a commercial microwave oven is discussed. This device when installed in conjunction with the existing circuitry of SHARP MICROWAVE OVEN (model R-9524) is capable of automatically advancing through a sequence of thawing recipes programmed and stored in the memory bank of the oven. The device therefore eliminates or minimizes human operator action needed in previous prototype version of a blood thawing device.

INTRODUCTION AND BACKGROUND

In today's society where progress in Medical Science has facilitated advancement in health care, there is increasing need for readily available biological tissues for transplantation. In the armed services, this need can be critical during times of confrontation and war. Presently, cryopreserved tissues such as blood are thawed upon need using a warm water bath. This method requires a water bath, heating element, pumps, filter and such other ancillary equipment which makes it cumbersome. Additionally, this may require transportation of clean water to far away territories to ensure safe, uncontaminated thawing. Over the last several years the U. S. Navy has supported research and development of an alternate method for thawing frozen blood components. Under that support this investigator had developed a working prototype to thaw frozen blood in bulk quantities using modified commercial microwave ovens. The evaluation tests conducted at NBRL-Boston demonstrated the success of this method. The thawed units although, were judged to be excellent from percent recovery and biochemical standpoints, concern was raised by the user (NBRL) over the lack of "user friendliness" of the system. Specifically, the operation of the system needed continued operation action and supervision. This was unacceptable for large-scale clinical usage intended by the Navy. Development of an electronic controller was therefore undertaken under the present contract. A control device was developed fabricated, tested and delivered to the Naval Blood Research

Laboratory (NEPL) Boston in September 1987 via the Naval Ocean Systems Center (NOSC) in San Diego.

DISCUSSION

The model R-9524 SHARP microwave oven is equipped with an auto-touch control panel shown in figure 1. The touch panel is hard wired to a proprietary micro processor (LSI-1x2029YA). The oven operation is controlled by touching appropriate pads arranged on the surface of the control panel. A lighted digital readout displays cooking time, temperature, memory sequence etc. during the operation of the oven.

As seen from figure 1, the memory bank has three memories (numbered Memory 1, 2 and 3). Each memory can be accessed independently and only through operator action. Each memory can be programmed to contain a maximum of three numerical settings corresponding to a temperature or time of cooking at a desired power level (high, medium high, medium, medium low or hold). Upon initiating the cooking sequence (started by touching COOK pad), the oven sequences through the numerical settings and performs the cooking operations within the memory. The oven stops after sequencing through a memory. Since the three memory banks are independent, the oven is not capable of accessing cooking instructions stored in subsequent memory banks without operator action. For example if the operator initiates cooking in Memory 1, the oven cannot access instructions in Memory 2 or Memory 3. Consequently, the oven stops at the conclusion of Memory 1.

This feature caused difficulty in thawing frozen blood. A safe blood thawing protocol calls for a recipe containing 7 to 9 numerical settings requiring all three memory banks. This meant that the use of Sharp Model R-9524 necessitated constant operator action to push COOK pad at the conclusion of instructions in each memory. To circumvent this difficulty an electronic controller was devised. This device facilitated accessing all three memories automatically without operator action.

THE AUTO CONTROLLER

The controller is schematized in fig. 2. It is designed to be operated in conjunction with the existing control circuitry of the oven. In principle it is a switch that activates the "next" memory upon completion of instructions in the "present" memory.

The device is activated by the audible signal given-out by the oven controller at the end of each memory sequence. The circuitry distinguishes this audio signal from other audio signals and produces a signal similar (in shape and size) to the one produced by pushing "COOK" pad. This signal then activates the relay controlling the power-on that depends on the activation of Memory 2 or Memory 3. Thus automatic accessing of the instructions contained in Memories 1 - 3 becomes possible.

The circuitry shown in fig. 2 was fabricated on a standard proto-board and installed in two microwave ovens supplied by NOSC. The successful operation of the oven using this controller was tested and demonstrated to NBRL personnel during September 1987. The ovens are presently in operation at NOSC, San Diego and NBRL Boston.

AUTO-TOUCH™ PANEL




Figure 1. Auto-Touch Panel of Sharp Microwave Oven Model R-9524.

The views and conclusions contained in this report are those of the contractors and should not be interpreted as representing the official policies, either expressed or implied, of the Naval Ocean Systems Center or the U.S. Government.

Approved for public release; distribution is unlimited.